

TakeOne®- the only aseptic sampling system of it's kind

06 October 2020 | Features

TakeOne® Aseptic sampling device is completely ready to use & have multiple, independent sampling lines

SARTORIUS

Simplifying Progress

Not all samples are equal – categorization of a sampling plan by assay type reduces costs, improves sampling results and optimizes capacity of the sampling solution. Samples are drawn throughout the entire Pharma | Biopharmaceutical process for different steps like compounding, holding stage, cell expansion, fermentation, buffer | media preparation, purification, final fill finish etc.

Features and benefits:

- Single use – eliminates cross contamination.
- Pre-assembled & ready to use – saving time and increasing safety.
- 100% leak tested and Gamma sterilized.
- Available tank fittings – 25mm InGold , 1.5" and 2" TC
- TakeOne® Flex – sampling from single use containers.
- Quick and easy disconnect.
- No sample overflow- spring return actuator.
- Application driven containers – Centrifuge tubes, PETG bottles, EVA & PE bags , needle free sampling etc.
- Manifold availability – expand sampling capacity.

Image not found or type unknown

